





Pantone
2716c

C	40
M	29
Y	0
K	0
R	159
G	174
B	229






PREVENTIVE CARE PROFILE PLUS

FOR THE QUANTITATIVE ANALYSIS OF:

albumin (ALB)	globulin (GLOB)*
alanine aminotransferase (ALT)	glucose (GLU)
alkaline phosphatase (ALP)	potassium (K+)
aspartate aminotransferase (AST)	sodium (NA+)
blood urea nitrogen (BUN)	total carbon dioxide (tCO2)
calcium (CA)	total bilirubin (TBIL)
chloride (CL-)	total protein (TP)
creatinine (CRE)	*calculated value



Caution


Sodium Heparin **MUST NOT** be used when collecting blood samples for use with this panel.

Bei der Entnahme von Blutproben für den Gebrauch mit diesem Panel **DARF KEIN** Natriumheparin verwendet werden.


NE PAS UTILISER l'héparine de sodium lors de la collecte des échantillons de sang qui seront utilisés avec ce bilan.

NO DEBE USARSE heparina sódica durante la recolección de muestras de sangre para ser utilizadas con este panel.


NON UTILIZZARE MAI eparina di sodio nella raccolta di campioni di sangue da utilizzare con questo pannello.




Abaxis, Inc.
3240 Whipple Road
Union City, CA 94587
800-822-2947




ABAXIS Europe GmbH
Bunsenstr. 9-11
64347 Griesheim
Germany
+49 6155 780 210



One Use




20 min



IVD

FOR VETERINARY USE ONLY

PN: 500-1047



2°C — 8°C

Made in USA

For patent information, see www.abaxis.com/about_us/patents.html

PN: 500-7183 Rev. A

